

# Agilis AG-UC2 & AG-UC8

## *Agilis™ Series Controllers*



**Newport®**

Experience | Solutions

## **Command Library API Manual**

**V2.0.0**

*For Motion, Think Newport™*

# Table of Contents

<b>1</b>	<b>INTRODUCTION.....</b>	<b>1</b>
1.1	PURPOSE .....	1
1.2	OVERVIEW .....	1
<b>2</b>	<b>COMMAND INTERFACE .....</b>	<b>1</b>
2.1	CONSTRUCTOR.....	1
2.2	FUNCTIONS .....	1
2.2.1	<i>General</i> .....	1
2.2.1.1	OpenInstrument.....	1
2.2.1.2	CloseInstrument.....	2
2.2.1.3	GetDevices .....	2
2.2.1.4	WriteToInstrument .....	2
2.2.2	<i>Commands</i> .....	2
2.2.2.1	CC_Get .....	2
2.2.2.2	CC_Set .....	3
2.2.2.3	JA_Get .....	3
2.2.2.4	JA_Set.....	3
2.2.2.5	MA.....	3
2.2.2.6	ML .....	4
2.2.2.7	MR.....	4
2.2.2.8	MV_Get .....	4
2.2.2.9	MV_Set.....	4
2.2.2.10	PA_Get.....	5
2.2.2.11	PA_Set.....	5
2.2.2.12	PH .....	5
2.2.2.13	PR.....	5
2.2.2.14	RS.....	6
2.2.2.15	ST .....	6
2.2.2.16	SU_Get.....	6
2.2.2.17	SU_Set .....	6
2.2.2.18	TE.....	7
2.2.2.19	TP.....	7
2.2.2.20	TS.....	7
2.2.2.21	VE.....	7
2.2.2.22	ZP .....	8
	<b>SERVICE FORM.....</b>	<b>9</b>

# Agilis AG-UC2 & AG-UC8

## Agilis™ Series Controllers

## 1 Introduction

---

### 1.1 Purpose

The purpose of this document is to describe the application programming interface (API) of the command library (AgilisCmdLib.dll) that is used to communicate with the AGILIS-UC2 or AGILIS-UC8 device.

### 1.2 Overview

The command library provides public methods to communicate with any AGILIS-UC2 or AGILIS-UC8 device and these methods work in both synchronous and asynchronous mode. Many of the ASCII commands that can be programmatically sent to the instrument have a corresponding method that can be called in the command library. For example, the ASCII “VE” command can be sent to the instrument to get the controller version, and the command library has a corresponding public method “VE” that returns the controller version and error information. For more information on a particular ASCII command see the manual for the controller.

## 2 Command Interface

---

### 2.1 Constructor

AgilisCmds ()

The constructor is used to create an instance of the command library.

### 2.2 Functions

#### 2.2.1 General

##### 2.2.1.1 OpenInstrument

###### Syntax

int OpenInstrument(string strDeviceKey)

string strDeviceKey: the device key is a serial COM port

return: 0 = successful or -1 = failure

###### Description

This function allows opening communication with the selected device. If the opening failed, the returned code is -1.

### 2.2.1.2 CloseInstrument

#### Syntax

int CloseInstrument()

return: 0 = successful or -1 = failure

#### Decription

This function allows closing communication with the selected device. If the closing failed, the returned code is -1.

### 2.2.1.3 GetDevices

#### Syntax

string[] GetDevices()

return: list of strings that contains the accessible COM ports.

#### Decription

This function returns the list of connected devices available to communicate.

### 2.2.1.4 WriteToInstrument

#### Syntax

int WriteToInstrument(string command, ref string resp, int stage)

command: Instrument command

resp: Response of the command

stage: Instrument Stage

return: function error

#### Decription

This overridden function Queries or writes the command given by the user to the instrument.

## 2.2.2 Commands

### 2.2.2.1 CC\_Get

#### Syntax

int CC\_Get(out int ChannelNumber, out string errstring)

ChannelNumber: ChannelNumber

errString: The failure reason

return: 0 in success and -1 on failure

#### Description

This function is used to process synchronous CC Get command which is used to Select channel (AGILIS-UC8 only). Refer to the AG-UC2-UC8 User's manual to get the command description.

### 2.2.2.2 CC\_Set

**Syntax**

int CC\_Set(int ChannelNumber, out string errstring)

ChannelNumber: ChannelNumber

errString: The failure reason

return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous CC Set command which is used to Select channel (AGILIS-UC8 only). Refer to the AG-UC2-UC8 User's manual to get the command description.

### 2.2.2.3 JA\_Get

**Syntax**

int JA\_Get(int controllerAddress, out int JogMode, out string errstring)

controllerAddress: controllerAddress

JogMode: JogMode

errString: The failure reason

return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous JA Get command which is used to Jog motion. Refer to the AG-UC2-UC8 User's manual to get the command description.

### 2.2.2.4 JA\_Set

**Syntax**

int JA\_Set(int controllerAddress, int JogMode, out string errstring)

controllerAddress: controllerAddress

JogMode: JogMode

errString: The failure reason

return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous JA Set command which is used to Jog motion. Refer to the AG-UC2-UC8 User's manual to get the command description.

### 2.2.2.5 MA

**Syntax**

int MA(int controllerAddress, out int Distance, out string errstring)

controllerAddress: controllerAddress

Distance: Distance

errString: The failure reason

return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous MA Get command which is used to Measure current position. Refer to the AG-UC2-UC8 User's manual to get the command description.

### 2.2.2.6 ML

#### **Syntax**

int ML(out string errstring)

errString: The failure reason

return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous ML Set command which is used to Local mode. Refer to the AG-UC2-UC8 User's manual to get the command description.

### 2.2.2.7 MR

#### **Syntax**

int MR(out string errstring)

errString: The failure reason

return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous MR Set command which is used to Remote mode. Refer to the AG-UC2-UC8 User's manual to get the command description.

### 2.2.2.8 MV\_Get

#### **Syntax**

int MV\_Get(int controllerAddress, out int JogMode, out string errstring)

controllerAddress: controllerAddress

JogMode: JogMode

errString: The failure reason

return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous MV Get command which is used to Move to limit. Refer to the AG-UC2-UC8 User's manual to get the command description.

### 2.2.2.9 MV\_Set

#### **Syntax**

int MV\_Set(int controllerAddress, int JogMode, out string errstring)

controllerAddress: controllerAddress

JogMode: JogMode

errString: The failure reason

return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous MV Set command which is used to Move to limit. Refer to the AG-UC2-UC8 User's manual to get the command description.

### 2.2.2.10 PA\_Get

**Syntax**

int PA\_Get(int controllerAddress, out int Target, out string errstring)

controllerAddress: controllerAddress

Target: Target

errString: The failure reason

return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous PA Get command which is used to Absolute move. Refer to the AG-UC2-UC8 User's manual to get the command description.

### 2.2.2.11 PA\_Set

**Syntax**

int PA\_Set(int controllerAddress, int Target, out string errstring)

controllerAddress: controllerAddress

Target: Target

errString: The failure reason

return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous PA Set command which is used to Absolute move. Refer to the AG-UC2-UC8 User's manual to get the command description.

### 2.2.2.12 PH

**Syntax**

int PH(out int LimitStatus, out string errstring)

LimitStatus: LimitStatus

errString: The failure reason

return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous PH Get command which is used to Tell limit status. Refer to the AG-UC2-UC8 User's manual to get the command description.

### 2.2.2.13 PR

**Syntax**

int PR(int controllerAddress, long NumberSteps, out string errstring)

controllerAddress: controllerAddress

NumberSteps: NumberSteps

errString: The failure reason

return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous PR Set command which is used to Relative move. Refer to the AG-UC2-UC8 User's manual to get the command description.

**2.2.2.14 RS****Syntax**

int RS(out string errstring)

errString: The failure reason

return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous RS Set command which is used to Reset. Refer to the AG-UC2-UC8 User's manual to get the command description.

**2.2.2.15 ST****Syntax**

int ST(int controllerAddress, out string errstring)

controllerAddress: controllerAddress

errString: The failure reason

return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous ST Set command which is used to Stop motion. Refer to the AG-UC2-UC8 User's manual to get the command description.

**2.2.2.16 SU\_Get****Syntax**

int SU\_Get(int controllerAddress, string Direction, out int StepAmplitude, out string errstring)

controllerAddress: controllerAddress

Direction: Direction

StepAmplitude: StepAmplitude

errString: The failure reason

return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous SU Get command which is used to Set step amplitude. Refer to the AG-UC2-UC8 User's manual to get the command description.

**2.2.2.17 SU\_Set****Syntax**

int SU\_Set(int controllerAddress, string Direction, int StepAmplitude, out string errstring)

controllerAddress: controllerAddress

Direction: Direction

StepAmplitude: StepAmplitude

errString: The failure reason



return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous SU Set command which is used to Set step amplitude. Refer to the AG-UC2-UC8 User's manual to get the command description.

### **2.2.2.18 TE**

#### **Syntax**

int TE(out int ErrorPreviousCommand, out string errstring)

ErrorPreviousCommand: ErrorPreviousCommand

errString: The failure reason

return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous TE Get command which is used to Get error of previous command. Refer to the AG-UC2-UC8 User's manual to get the command description.

### **2.2.2.19 TP**

#### **Syntax**

int TP(int controllerAddress, out int NumberSteps, out string errstring)

controllerAddress: controllerAddress

NumberSteps: NumberSteps

errString: The failure reason

return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous TP Get command which is used to Tell number of steps. Refer to the AG-UC2-UC8 User's manual to get the command description.

### **2.2.2.20 TS**

#### **Syntax**

int TS(int controllerAddress, out int AxisStatus, out string errstring)

controllerAddress: controllerAddress

AxisStatus: AxisStatus

errString: The failure reason

return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous TS Get command which is used to Get axis status. Refer to the AG-UC2-UC8 User's manual to get the command description.

### **2.2.2.21 VE**

#### **Syntax**

int VE(out string ControllerVersion, out string errstring)

ControllerVersion: ControllerVersion

errString: The failure reason

return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous VE Get command which is used to Firmware version. Refer to the AG-UC2-UC8 User's manual to get the command description.

### **2.2.2.22 ZP**

#### **Syntax**

int ZP(int controllerAddress, out string errstring)

controllerAddress: controllerAddress

errString: The failure reason

return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous ZP Set command which is used to Zero position. Refer to the AG-UC2-UC8 User's manual to get the command description.

# Service Form

## Your Local Representative

Tel.: \_\_\_\_\_

Fax: \_\_\_\_\_

Name: \_\_\_\_\_

Company: \_\_\_\_\_

Address: \_\_\_\_\_

Country: \_\_\_\_\_

P.O. Number: \_\_\_\_\_

Item(s) Being Returned: \_\_\_\_\_

Model#: \_\_\_\_\_

Return authorization #: \_\_\_\_\_

*(Please obtain prior to return of item)*

Date: \_\_\_\_\_

Phone Number: \_\_\_\_\_


Fax Number: \_\_\_\_\_

Serial #: \_\_\_\_\_

Description: \_\_\_\_\_

Reasons of return of goods (please list any specific problems): \_\_\_\_\_

[illegible]



**Newport®**

Experience | Solutions

Visit Newport Online at:  
**[www.newport.com](http://www.newport.com)**

**North America & Asia**

Newport Corporation  
1791 Deere Ave.  
Irvine, CA 92606, USA

**Sales**

Tel.: (800) 222-6440  
e-mail: [sales@newport.com](mailto:sales@newport.com)

**Technical Support**

Tel.: (800) 222-6440  
e-mail: [tech@newport.com](mailto:tech@newport.com)

**Service, RMAs & Returns**

Tel.: (800) 222-6440  
e-mail: [rma.service@newport.com](mailto:rma.service@newport.com)

**Europe**

MICRO-CONTROLE Spectra-Physics S.A.S  
1, rue Jules Guesde – Bât. B  
ZI Bois de l'Épine – BP189  
91006 Evry Cedex  
France

**Sales**

Tel.: +33 (0)1.60.91.68.68  
e-mail: [france@newport-fr.com](mailto:france@newport-fr.com)

**Technical Support**

e-mail: [tech\\_europe@newport.com](mailto:tech_europe@newport.com)

**Service & Returns**

Tel.: +33 (0)2.38.40.51.55